

**SP-2 拉酷**

**LUCK**

2-Layer Composite Resin Teeth

 Certified by  
SGS  
ISO13485:2003



# LUCK SP-2 RESIN TEETH

## ▼ TAPERING UPPER

<b>T11</b> A: 11.9 B: 9.6 C: 51.6	
<b>T12</b> A: 11.3 B: 9.2 C: 49.1	
<b>T13</b> A: 10.8 B: 8.7 C: 46.6	
<b>T14</b> A: 10.2 B: 8.2 C: 44.2	
<b>T15</b> A: 9.8 B: 7.8 C: 42.2	

## 混合形 COMBINATION UPPER

<b>R5</b> A: 11.5 B: 9.1 C: 50.5	
<b>R4</b> A: 10.9 B: 8.7 C: 48.1	
<b>R3</b> A: 10.7 B: 8.5 C: 47.2	
<b>R2</b> A: 10.4 B: 8.2 C: 45.7	
<b>R1</b> A: 10.0 B: 7.9 C: 43.9	

## ■ SQUARE UPPER

<b>B1</b> A: 11.0 B: 9.2 C: 49.7	
<b>B2</b> A: 10.4 B: 8.7 C: 47.2	
<b>B3</b> A: 9.9 B: 8.3 C: 44.7	
<b>B4</b> A: 9.3 B: 7.8 C: 42.3	
<b>B5</b> A: 8.8 B: 7.3 C: 39.7	

## ■ SHORT SQUARE UPPER


<b>V5</b> A: 11.4 B: 9.1 C: 48.9	
<b>V4</b> A: 10.8 B: 8.6 C: 46.1	
<b>V3</b> A: 10.2 B: 8.1 C: 43.8	
<b>V2</b> A: 9.7 B: 7.7 C: 41.6	
<b>V1</b> A: 9.1 B: 7.3 C: 39.2	

# LUCK SP-2 RESIN TEETH

COMBINATION	UPPER
<b>TO5</b> A: 12.8 B: 9.0 C: 46.7	
<b>TO4</b> A: 12.1 B: 8.5 C: 44.2	
<b>TO3</b> A: 11.5 B: 7.9 C: 40.8	
<b>TO2</b> A: 10.9 B: 7.5 C: 39.0	
<b>TO1</b> A: 10.5 B: 7.2 C: 37.4	

SQUARE	LOWER
<b>L5</b> A: 10.5 B: 5.5 C: 37.5	
<b>L4</b> A: 9.6 B: 5.3 C: 36.2	
<b>L3</b> A: 9.2 B: 5.0 C: 33.6	
<b>L2</b> A: 8.7 B: 4.7 C: 31.8	
<b>L1</b> A: 8.2 B: 4.5 C: 30.8	

TAPERING	LOWER
<b>L11</b> A: 11.8 B: 5.3 C: 36.2	
<b>L12</b> A: 11.2 B: 5.1 C: 34.4	
<b>L13</b> A: 10.6 B: 4.8 C: 32.6	
<b>L14</b> A: 10.1 B: 4.5 C: 30.8	
<b>L15</b> A: 9.5 B: 4.3 C: 29.0	



## Reading Guide


- A: Length of central cutting teeth (excluding neck)
- B: Width of central cutting teeth
- C: Width of 6 anterior teeth

The scale of size in chart to real is 1:1  
All data is measured in millimetre (mm)


# LUCK SP-2 RESIN TEETH

POSTERIOR MOLD (18°) for BUCCAL Occlusion


M5




M4



M3



M2



2-Layer Acrylic Teeth Combination Table

Upper	Lower	Posterior	Upper	Lower	Posterior	Upper	Lower	Posterior
T11	L11	M5/M34	V1	L1	M2/M28	TO1	L15	M2/M28
T12	L12	M4/M32	V2	L2	M3/M30	TO2	L14	M3/M30
T13	L13	M4/M32	V3	L3	M3/M30	TO3	L13	M3/M30
T14	L14	M3/M30	V4	L4	M4/M32	TO4	L12	M4/M32
T15	L15	M2/M28	V5	L5	M4/M32	TO5	L11	M4/M32
B1	L5	M4/M32	R1	L15	M2/M28			
B2	L4	M4/M32	R2	L14	M3/M30			
B3	L3	M3/M30	R3	L13	M3/M30			
B4	L2	M3/M30	R4	L12	M4/M32			
B5	L1	M2/M28	R5	L11	M4/M32			


Remark:

Combinations offered here are for reference only. User can choose their combinations according to clinical needs. Combinations provided here are subject to changes on our website.


# LUCK SP-2 RESIN TEETH

POSTERIOR MOLD (28°) for LINGUAL Occlusion


M34




M32



M30



M28



**SHANGHAI PIGENG DENTAL MFG. CO., LTD.**

ADD. No.233 Fengyong Road, Qingcun, Fengxian, Shanghai, China 201414

TEL: 021-68880119 FAX: 021-68880117

E-mail:pigeon@pigeondental.com

Agent.