



PIGEON DENTAL
Make Confidence

SP-3 莎美特

SUMMIT

3-Layer Composite Resin Teeth

SP-3 玛露姬

MARUJI

3-Layer Acrylic Teeth



Cerificated by
SGS
ISO13485:2003



JIS T6506

SP-3 SUMMIT & MARUJI

▼ 尖形 TAPERING

UPPER

T5

A: 11.5
B: 9.0
C: 46.9



T4

A: 11.0
B: 8.5
C: 44.6



T3

A: 10.5
B: 8.0
C: 42.6



T2

A: 10.5
B: 7.5
C: 39.3



T1

A: 9.2
B: 7.2
C: 38.0



● 圓形 OVOID

UPPER

O5

A: 11.0
B: 9.0
C: 48.5



O4

A: 10.5
B: 8.5
C: 45.5



O3

A: 10.2
B: 8.0
C: 42.9



O2

A: 9.6
B: 7.6
C: 40.7



O1

A: 9.1
B: 7.2
C: 38.7



■ 方形 SQUARE

UPPER

S5

A: 11.5
B: 9.0
C: 49.5



S4

A: 11.0
B: 8.5
C: 47.5



S3

A: 10.5
B: 8.1
C: 44.7



S2

A: 10.0
B: 7.6
C: 42.4



S1

A: 9.5
B: 7.3
C: 40.0



■ 短方形 SHORT SQUARE

UPPER

SS5

A: 11.7
B: 9.0
C: 47.6



SS4

A: 10.2
B: 8.4
C: 44.6



SS3

A: 10.2
B: 8.0
C: 42.2



SS2

A: 9.6
B: 7.5
C: 39.6



SS1

A: 9.1
B: 7.2
C: 37.6



SP-3 SUMMIT & MARUJI



全口牙齿组合表 3-Layer Acrylic Teeth Combination Table

Upper	Lower	Posterier	Upper	Lower	Posterier	Upper	Lower	Posterier	Upper	Lower	Posterier
T5	LT5	M4/M34	O5	L5	M4/M34	S5	L5	M4/M34	SS5	L5	M4/M34
T4	LT4	M4/M32	O4	L4	M4/M32	S4	L4	M4/M33	SS4	L4	M4/M32
T3	LT3	M3/M30	O3	L3	M3/M30	S3	L3	M3/M30	SS3	L3	M3/M30
T2	LT2	M3/M30	O2	L2	M3/M30	S2	L2	M3/M30	SS2	L2	M2/M30
T1	LT1	M2/M28	O1	L1	M2/M28	S1	L1	M2/M28	SS1	L1	M2/M28


以上组合仅供参考，用户可以根据临床情况另行配比！

Remark:

Combinations offered here are for reference only. User can choose their combinations according to clinical needs.

Combinations provided here are subject to changes on our website.

尺寸单位：mm
图片牙齿比例：1:1







前牙尺寸：

- A - 中切牙长度
- B - 中切牙宽度
- C - 前牙6颗牙总宽度

Reading Guide

- A: Length of central cutting teeth (excluding neck)
- B: Width of central cutting teeth
- C: Width of 6 anterior teeth

The scale of size in chart to real is 1:1
All data is measured in millimetre (mm)



磨牙的排牙技巧：


- 1、舌侧集中咬合排牙法适用于欧美人群
- 2、颊侧集中咬合排牙法适用于亚洲人群
- 3、对于牙槽骨吸收较大的患者来说，选择颊侧集中排牙法最为理想！

SP-3 SUMMIT & MARUJI


POSTERIOR MOLD (28°) for LINGUALIZED Occlusion

此款磨牙适合舌集中排牙法


M34




M32



M30



M28




SP-3 SUMMIT & MARUJI


POSTERIOR MOLD (33°) for FULL BALANCE Occlusion

此款磨牙适合全平衡排牙法


M34




M32



M30



M28



上海贝琮齿材股份有限公司

SHANGHAI PIGENG DENTAL MFG. CO., LTD.

ADD. No.233 Fengyong Road, Qingcun, Fengxian,
Shanghai, China 201414

TEL: 021-68880119 FAX: 021-68880117

E-mail:pigeon@pigeondental.com

Agent.